
無創性血糖測量儀 NON-INVASIVE BLOOD GLUCOSE METER

利用近紅外線技術，以非創傷性方法測量糖尿病患者之血糖水平，從而鼓勵患者密切監察自己的血糖水平。
By using NIR to measure blood glucose non-invasively that encourages diabetic patients closely monitor their blood glucose levels.

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最近在臨床實驗中使用六個月特定的波長及交互確效校正效果
Calibrated results using an "ideal" pre-classifier.
The Pearson correlation coefficient is 0.86 (with $p < 0.01$).

Patients with Diabetes Mellitus are recommended to monitor their blood glucose levels frequently and regularly. Conventional blood glucose level tests involve taking blood samples which are difficult and painful to perform. Non-invasive Blood Glucose Meter applied the NIR technology, which measures blood glucose level. It helps reducing sufferings, the size and cost of the corresponding hardware devices to be developed, and facilitates portable and continuous monitoring of blood glucose level. The measurements could be digitized and stored as health records, and such continuous monitoring health records could then be retrieved easily for health advices.

Principal Investigator
Prof. Joanne Wai-yea Chung
School of Nursing
Contact Details

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Patent Application No: PCT/IB2006/001142(PCT) 11/122,325(USA)

Special Features and Advantages

The new method is needleless and painless, neither implanting a sensor subcutaneously nor pricking fingers for a blood sample. Thus, there is zero risk of infection. It is cost saving because no laboratory tests or consumables, e.g. testing reagent are needed. The measuring method is extremely user-friendly, which is to put the thumb over the NIR sensor for 5 to 6 seconds.

By identifying and reducing the number of wavelengths required, this method can reduce the size of a thimble.

Application

The method not only facilitates blood glucose measurements, it may also be extended to the monitoring of other blood substrates such as cholesterol and lactic acid.

Awards

- Gold Medal with the Congratulations of the Jury, 35th International Exhibition of Inventions New Techniques & Products, Geneva (April 2007)
- Bronze Award, Brussels Eureka 2005 -the 54th World Exhibition of Innovation, Research and New Technology (November 2005)